

Disposable Spo2 Sensors



Features

- *Engineered to function equivalent to OEMs*
- *Quality approved with CE Certificate and FDA 510K*
- *Single patient use eliminates cross-contamination issues*
- *Reliable performance with sensitive and accurate responses*
- *Cost effective*

Available across the entire line of primary brand compatible disposable sensors

- | | | | |
|------------------|--------------------|--------------|-----------------------|
| ▪ Nellcor | ▪ Masimo | ▪ GE | ▪ Nihon Kohden |
| ▪ Nonin | ▪ BCI/Smith | ▪ CSI | |



PN: SA-001X-XC
Price: \$20.00/pc

Oximax Adapter Cable


The Oximax Adapter Cable cuts the cost of Oximax disposable sensors allowing you to use the non-Oximax Nellcor disposable sensors with same performance functionally as the Oximax sensors.

Nellcor- Oximax


Cable length: 0.9m

Package: 24 pcs/box

Price: USD192.00/box



SSD-001XC-W09A



SSD-001XC-S09PI



Sensor Connector


| Part Number | SSD-001XC-W09A | SSD-001XC-W09P | SSD-001XC-W09I | SSD-001XC-W09AN | SSD-001XC-S09AN | SSD-001XC-S09PI |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM# | MAXA, MAXAL | MAXP | MAXI | MAXAN | MAXA, MAXI | MAXP, MAXI |
| Patient Sizes | Adult >40kg | Pediatric 10-50kg | Infant 3-20kg | Neonate/Adult <3kg or >40kg | Adult/Neonate >40kg or <3kg | Pediatric/Infant 3-40kg |
| Tapes | Adhesive Non-woven fiber | | | | Non-adhesive foam | |

Nellcor- Oxismart (non Oximax)


Cable length: 0.9m

Package: 24 pcs/box


Price: USD168.00/box



SSD-001-W09A



SSD-001-S09AN



Sensor Connector


| Part Number | SSD-001-W09A | SSD-001-W09P | SSD-001-W09I | SSD-001-W09AN | SSD-001-S09AN | SSD-001-S09PI |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM# | D25 | D20 | I20 | N25 | N25, D25 | D20, I20 |
| Patient Sizes | Adult >40kg | Pediatric 10-50kg | Infant 3-20kg | Neonate/Adult <3kg or >40kg | Adult/Neonate >40kg or <3kg | Pediatric/Infant 3-40kg |
| Tapes | Adhesive Non-woven fiber | | | | Non-adhesive foam | |

Masimo- LNCS

Cable length: 0.9m

Package: 20 pcs/box



Price: USD180.00/box



SSD-005RB-W09A



SSD-005RB-S09AN



Sensor Connector

| Part Number | SSD-005RB-W09A | SSD-005RB-W09P | SSD-005RB-W09I | SSD-005RB-W09AN | SSD-005RB-S09AN | SSD-005RB-S09PI |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM# | 1859,2317 | 1860,2318 | 1861,2328 | 2329,1862 | 2606 | 1860,2328 |
| Patient Sizes | Adult >40kg | Pediatric 10-50kg | Infant 3-20kg | Neonate/Adult <3kg or >40kg | Adult/Neonate >40kg or <3kg | Pediatric/Infant 3-40kg |
| Tapes | Adhesive Non-woven fiber | | | | Non-adhesive foam | |

Masimo- M-LNCS

Cable length: 0.9m

Package: 20 pcs/box



Price: USD220.00/box



SSD-007RB-W09A



SSD-007RB-S09AN



Sensor Connector


| Part Number | SSD-007RB-W09A | SSD-007RB-W09P | SSD-007RB-W09I | SSD-007RB-W09AN | SSD-007RB-S09AN | SSD-007RB-S09PI |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM# | 2805,2509 | 2510,2511 | 2512,2513 | 2514,2515 | 2805,2509 | 2510,2511 |
| Patient Sizes | Adult >40kg | Pediatric 10-50kg | Infant 3-20kg | Neonate/Adult <3kg or >40kg | Adult/Neonate >40kg or <3kg | Pediatric/Infant 3-40kg |
| Tapes | Adhesive Non-woven fiber | | | | Non-adhesive foam | |

Masimo- LNOP

Cable length: 0.9m

Package: 20 pcs/box


Price: USD200.00/box



SSD-006RB-W09A



SSD-006RB-S09AN



Sensor Connector


| Part Number | SSD-006RB-W09A | SSD-006RB-W09P | SSD-006RB-W09I | SSD-006RB-W09AN | SSD-006RB-S09AN | SSD-006RB-S09PI |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM# | 1001 | 1025 | 1800 | 2358 | 1001 | 1025 |
| Patient Sizes | Adult >40kg | Pediatric 10-50kg | Infant 3-20kg | Neonate/Adult <3kg or >40kg | Adult/Neonate >40kg or <3kg | Pediatric/Infant 3-40kg |
| Tapes | Adhesive Non-woven fiber | | | | Non-adhesive foam | |

GE- Datex

Cable length: 0.9m

Package: 24 pcs/box


Price: USD180.00/box



SSD-013-W09A



SSD-013-S09AN



Sensor Connector


| Part Number | SSD-013-W09A | SSD-013 -W09P | SSD-013-W09I | SSD-013-W09AN | SSD-013-S09AN | SSD-013-S09PI |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM# | SAS-AF | SAS-AP | SAS-AP | SAS-AF | SAS-AF | SAS-AP |
| Patient Sizes | Adult >40kg | Pediatric 10-50kg | Infant 3-20kg | Neonate/Adult <3kg or >40kg | Adult/Neonate >40kg or <3kg | Pediatric/Infant 3-40kg |
| Tapes | Adhesive Non-woven fiber | | | | Non-adhesive foam | |

GE- TruSignal

Cable length: 0.9m

Package: 24 pcs/box


Price: USD200.00/box



SSD-018-W09A



SSD-018-S09AN



Sensor Connector


| Part Number | SSD-018-W09A | SSD-018-W09P | SSD-018-W09I | SSD-018-W09AN | SSD-018-S09AN | SSD-018-S09PI |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM# | TS-AF, TS-AAW | TS-AP, TS-PAW | TS-AF | TS-AF | TS-AF, TS-AAW | TS-AP, TS-PAW |
| Patient Sizes | Adult >40kg | Pediatric 10-50kg | Infant 3-20kg | Neonate/Adult <3kg or >40kg | Adult/Neonate >40kg or <3kg | Pediatric/Infant 3-40kg |
| Tapes | Adhesive Non-woven fiber | | | | Non-adhesive foam | |

GE- OxyTip

Cable length: 0.9m

Package: 24 pcs/box


Price: USD200.00/box



SSD-019-W09A



SSD-019-S09AN



Sensor Connector


| Part Number | SSD-019-W09A | SSD-019-W09P | SSD-019-W09I | SSD-019-W09AN | SSD-019-S09AN | SSD-019-S09PI |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM# | OXY-AF | OXY-AP,OXY-AF | OXY-AF | OXY-AF | OXY-AF | OXY-AP,OXY-AF |
| Patient Sizes | Adult >40kg | Pediatric 10-50kg | Infant 3-20kg | Neonate/Adult <3kg or >40kg | Adult/Neonate >40kg or <3kg | Pediatric/Infant 3-40kg |
| Tapes | Adhesive Non-woven fiber | | | | Non-adhesive foam | |

Nihon Kohden


Cable length: 0.9m

Package: 24 pcs/box


Price: USD180.00/box



SSD-043-W09A



SSD-043-S09AN



Sensor Connector


| Part Number | SSD-043-W09A | SSD-043-W09P | SSD-043-W09I | SSD-043-W09AN | SSD-043-S09AN | SSD-043-S09PI |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM# | TL-271 | TL-272 | TL-274 | TL-273 | TL-273 | TL-272/T274 |
| Patient Sizes | Adult >40kg | Pediatric 10-50kg | Infant 3-20kg | Neonate/Adult <3kg or >40kg | Adult/Neonate >40kg or <3kg | Pediatric/Infant 3-40kg |
| Tapes | Adhesive Non-woven fiber | | | | Non-adhesive foam | |

Nonin

Cable length: 0.9m

Package: 24 pcs/box


Price: USD168.00/box



SSD-040-W09A



SSD-040-S09PI



Sensor Connector


| Part Number | SSD-040-W09A | SSD-040-W09P | SSD-040-W09I | SSD-040-W09AN | SSD-040-S09AN | SSD-040-S09PI |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM# | 7000A, 6000CA | 7000P, 6000CP | 7000I, 6000CI | 7000N, 6000CN | 7000A, 7000N, | 7000P, 7000I, |
| Patient Sizes | Adult >40kg | Pediatric 10-50kg | Infant 3-20kg | Neonate/Adult <3kg or >40kg | Adult/Neonate >40kg or <3kg | Pediatric/Infant 3-40kg |
| Tapes | Adhesive Non-woven fiber | | | | Non-adhesive foam | |

BCI/Smith

Cable length: 0.9m

Package: 24 pcs/box



Price: USD168.00/box



SSD-036-W09A



SSD-036-S09AN



Sensor Connector

| Part Number | SSD-036-W09A | SSD-036-W09P | SSD-036-W09I | SSD-036-W09AN | SSD-036-S09AN | SSD-036-S09PI |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM# | 1300 | 1301 | 1303 | 1302 | 1300, 1302 | 1301, 1303 |
| Patient Sizes | Adult >40kg | Pediatric 10-50kg | Infant 3-20kg | Neonate/Adult <3kg or >40kg | Adult/Neonate >40kg or <3kg | Pediatric/Infant 3-40kg |
| Tapes | Adhesive Non-woven fiber | | | | Non-adhesive foam | |

CSI (Criticare)

Cable length: 0.9m

Package: 24 pcs/box

Price: USD180.00/box



Sensor Connector

| Part Number | SSD-039-W09A | SSD-039-W09P | SSD-039-W09I | SSD-039-W09AN | SSD-039-S09AN | SSD-039-S09PI |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM# | 570SD | 571SD | 572SD | 573SD | 570SD | 571SD, 572SD |
| Patient Sizes | Adult >40kg | Pediatric 10-50kg | Infant 3-20kg | Neonate/Adult <3kg or >40kg | Adult/Neonate >40kg or <3kg | Pediatric/Infant 3-40kg |
| Tapes | Adhesive Non-woven fiber | | | | Non-adhesive foam | |