

Disposable Temperature Probes



Features

- *Engineered to function equivalent to OEMs*
- *Quality approved with CE Certificate and FDA 510K*
- *Single patient use eliminates cross-contamination issues*
- *Reliable performance with sensitive and accurate responses*
- *Cost effective*



Available across the entire line of primary brand compatible disposable sensors

▪ **Philips**

▪ **GE**

▪ **YSI400**

Philips- 400 series

Disposable temperature probes


Cable length: 0.9m

Package: 25 pcs/box

Price: USD125.00/box



TSD-PM-GS09



TSD-PM-AG09



Probe Connector

Part Number	Descriptions	OEM Part number
TSD-PM-GS09	Skin temperature probe, disposable	21091A,989803100951
TSD-PM-AG09	Adult general temperature probe, disposable	989803105321 21090A,989803100941

GE- 400 series

Disposable temperature probes

Cable length: 0.9m

Package: 25 pcs/box


Price: USD125.00/box



TSD-YSI-GS09



TSD-YSI-AG09



Probe Connector

Part Number	Descriptions	OEM Part number
TSD-YSI-GS09	Skin temperature probe, disposable	M1024222 (<i>sterile</i>) 1017EU (<i>non-sterile</i>)
TSD-YSI-AG09	Adult general temperature probe, disposable	M1024229 (<i>sterile,9t</i>) M1024231 (<i>sterile,12f</i>)


Adapter cables for disposable probes

Cable length: 3.0m

Price: USD




TA-PH-PM30



TA-MRA-YSI30



used with philips compatible disposable probes



used with GE compatible disposable probes

Manufacturer	Monitor End	Probe end	Part Number	OEM Reference	Price
YSI			TA-Y400-PM30		\$25.00
			TA-Y400-YSI30	165640(1.3m) 165641(2.8m)	\$25.00
GE			TA-MQ-PM30		\$25.00
			TA-MQ-YSI30	2021701-001(3.6m) 2021702-002(1.2m)	\$25.00
			TA-MQ-2YSI30	2021704-002(3.6m) 2021705-002(1.2m)	\$25.00
Infinium			TA-IFN-PM30		\$25.00
			TA-IFN-YSI30		\$25.00
Mindray			TA-MRA-PM30		\$25.00
			TA-MRA-YSI30	040-000091-00	\$25.00
Mindray-Datascope			TA-MRB-PM30		\$25.00
			TA-MRB-YSI30	040-001235-00	\$25.00
Philips			TA-PH-PM30	21082A/989803100921 21082B/989803100931	\$25.00
			TA-PH-YSI30		\$30.00
Spacelabs			TA-SL-PM30		\$30.00
			TA-SL-YSI30		\$30.00