

# Single Patient Use ECG Leadwires



## Features

- *Quality approved with CE mark and FDA 510K*
- *Single Patient Use, eliminates cross-contamination issues*
- *Directly used with telemetry or OEM monitor trunk cables*
- *Cost effective and time saving*

Wide range selection of compatible Single Patient Use leadwires for the primary Brands:

- |                       |                                 |                     |
|-----------------------|---------------------------------|---------------------|
| ▪ <b>GE</b>           | ▪ <b>Kendall/CardinalHealth</b> | ▪ <b>Mindray</b>    |
| ▪ <b>Nihon Kohden</b> | ▪ <b>Philips</b>                | ▪ <b>Safety Din</b> |



## Conversion Cables and Adapters

Conversion cables and adapters are designed to connect into an existing monitor and change the "yoke" end to a common connection.

## Philips

Part Number	Product Descriptions	OEM Cross Reference	Packaging	Price
<b>1</b> ELD-PH3-90S-A	<b>IntelliVue</b> 3 lead, IntelliVue, Snap, 0.9m, AHA		10 pcs/box	
ELD-PH3-90S-I	<b>IntelliVue</b> 3 lead, IntelliVue, Snap, 0.9m, IEC		10 pcs/box	
ELD-PH3-90G-A	<b>IntelliVue</b> 3 lead, IntelliVue, Grabber, 0.9m, AHA	989803173121	10 pcs/box	
ELD-PH3-90G-I	<b>IntelliVue</b> 3 lead, IntelliVue, Grabber, 0.9m, IEC	989803174201	10 pcs/box	
ELD-PH5-90S-A	<b>IntelliVue</b> 5 lead, IntelliVue, Snap, 0.9m, AHA		10 pcs/box	
ELD-PH5-90S-I	<b>IntelliVue</b> 5 lead, IntelliVue, Snap, 0.9m, IEC		10 pcs/box	
<b>2</b> ELD-PH5-90G-A	<b>IntelliVue</b> 5 lead, IntelliVue, Grabber, 0.9m, AHA	989803173131	10 pcs/box	
<b>3</b> ELD-PH5-90G-I	<b>IntelliVue</b> 5 lead, IntelliVue, Grabber, 0.9m, IEC	989803174211	10 pcs/box	
ELD-PHV5-90S-A	<b>IntelliVue</b> 5 lead, IntelliVue, Snap, 0.9m, AHA, Chest lead		10 pcs/box	
ELD-PHV5-90S-I	<b>IntelliVue</b> 5 lead, IntelliVue, Snap, 0.9m, IEC, Chest lead		10 pcs/box	
ELD-PHV5-90G-A	<b>IntelliVue</b> 5 lead, IntelliVue, Grabber, 0.9m, AHA, Chest lead	989803192141	10 pcs/box	
ELD-PHV5-90G-I	<b>IntelliVue</b> 5 lead, IntelliVue, Grabber, 0.9m, IEC, Chest lead	989803192151	10 pcs/box	



## Philips

Part Number	Product Descriptions	OEM Cross Reference	Packaging	Price
ELD-TAA3-90S-A	<b>M2601A/M2600A/M2636A</b> 3 lead, Twin Pin, Snap, 0.9m, AHA		10 pcs/box	
ELD-TAA3-90S-I	<b>M2601A/M2600A/M2637A</b> 3 lead, Twin Pin, Snap, 0.9m, IEC		10 pcs/box	
ELD-TAA3-90G-A	<b>M2601A/M2600A/M2638A</b> 3 lead, Twin Pin, Grabber, 0.9m, AHA		10 pcs/box	
ELD-TAA3-90G-I	<b>M2601A/M2600A/M2639A</b> 3 lead, Twin Pin, Grabber, 0.9m, IEC		10 pcs/box	
ELD-TAA5-90S-A	<b>M2601A/M2600A/M2640A</b> 5 lead, Twin Pin, Snap, 0.9m, AHA	33113	10 pcs/box	
ELD-TAA5-90S-I	<b>M2601A/M2600A/M2641A</b> 5 lead, Twin Pin, Snap, 0.9m, IEC	33113E	10 pcs/box	
<b>1</b> ELD-TAA5-90G-A	<b>M2601A/M2600A/M2642A</b> 5 lead, Twin Pin, Grabber, 0.9m, AHA		10 pcs/box	
ELD-TAA5-90G-I	<b>M2601A/M2600A/M2643A</b> 5 lead, Twin Pin, Grabber, 0.9m, IEC		10 pcs/box	
<b>2</b> ELD-TPH3-90S-A	<b>M2601B/M4841A</b> 3 lead, IntelliVue, Snap, 0.9m, AHA		10 pcs/box	
ELD-TPH3-90S-I	<b>M2601B/M4842A</b> 3 lead, IntelliVue, Snap, 0.9m, IEC		10 pcs/box	
ELD-TPH3-90G-A	<b>M2601B/M4843A</b> 3 lead, IntelliVue, Grabber, 0.9m, AHA	989803173141	10 pcs/box	
ELD-TPH3-90G-I	<b>M2601B/M4844A</b> 3 lead, IntelliVue, Grabber, 0.9m, IEC		10 pcs/box	
ELD-TPH5-90S-A	<b>M2601B/M4845A</b> 5 lead, IntelliVue, Snap, 0.9m, AHA	33114	10 pcs/box	
ELD-TPH5-90S-I	<b>M2601B/M4846A</b> 5 lead, IntelliVue, Snap, 0.9m, IEC		10 pcs/box	
<b>3</b> ELD-TPH5-90G-A	<b>M2601B/M4847A</b> 5 lead, IntelliVue, Grabber, 0.9m, AHA	989803173151	10 pcs/box	
ELD-TPH5-90G-I	<b>M2601B/M4848A</b> 5 lead, IntelliVue, Grabber, 0.9m, IEC		10 pcs/box	



## Philips

Part Number	Product Descriptions	OEM Cross Reference	Packaging	Price
ELD-TMX3-90S-A	<b>MX40</b> 3 lead, Snap, 0.9m, AHA		10 pcs/box	
ELD-TMX3-90S-I	<b>MX40</b> 3 lead, Snap, 0.9m, IEC		10 pcs/box	
ELD-TMX3-90G-A	<b>MX40</b> 3 lead, Grabber, 0.9m, AHA		10 pcs/box	
ELD-TMX3-90G-I	<b>MX40</b> 3 lead, Grabber, 0.9m, IEC		10 pcs/box	
ELD-TMX5-90S-A	<b>MX40</b> 5 lead, Snap, 0.9m, AHA	33147	10 pcs/box	
ELD-TMX5-90S-I	<b>MX40</b> 5 lead, Snap, 0.9m, IEC		10 pcs/box	
<b>1</b> ELD-TMX5-90G-A	<b>MX40</b> 5 lead, Grabber, 0.9m, AHA	989803172031	10 pcs/box	
ELD-TMX5-90G-I	<b>MX40</b> 5 lead, Grabber, 0.9m, IEC		10 pcs/box	
ELD-TMX5S09AN-90S-A	<b>MX40</b> 5 lead, Snap, 0.9m, with Spo2 sensor, AHA	33149	10 pcs/box	
ELD-TMX5S09AN-90S-I	<b>MX40</b> 5 lead, Snap, 0.9m, with Spo2 sensor, IEC		10 pcs/box	
<b>2</b> ELD-TMX5S09AN-90G-A	<b>MX40</b> 5 lead, Grabber, 0.9m, with Spo2 sensors AHA	989803172051	10 pcs/box	
ELD-TMX5S09AN-90G-I	<b>MX40</b> 5 lead, Grabber, 0.9m, with Spo2 sensors IEC	989803172151	10 pcs/box	



## GE Healthcare

Part Number	Product Descriptions	OEM Cross Reference	Packaging	Price
ELD-MQ3-90S-A	<b>Multi-Link</b> 3 lead, Multi-Link, Snap, 0.9m, AHA		10 pcs/box	
ELD-MQ3-90S-I	<b>Multi-Link</b> 3 lead, Multi-Link, Snap, 0.9m, IEC	2052104-001	10 pcs/box	
<b>1</b> ELD-MQ3-90G-A	<b>Multi-Link</b> 3 lead, Multi-Link, Grabber, 0.9m, AHA		10 pcs/box	
ELD-MQ3-90G-I	<b>Multi-Link</b> 3 lead, Multi-Link, Grabber, 0.9m, IEC	2052133-001	10 pcs/box	
ELD-MQ5-90S-A	<b>Multi-Link</b> 5 lead, Multi-Link, Snap, 0.9m, AHA	33112	10 pcs/box	
ELD-MQ5-90S-I	<b>Multi-Link</b> 5 lead, Multi-Link, Snap, 0.9m, IEC	2052104-002	10 pcs/box	
<b>2</b> ELD-MQ5-90G-A	<b>Multi-Link</b> 5 lead, Multi-Link, Grabber, 0.9m, AHA		10 pcs/box	
ELD-MQ5-90G-I	<b>Multi-Link</b> 5 lead, Multi-Link, Grabber, 0.9m, IEC	2052133-002	10 pcs/box	
ELD-MQV5-90S-A	<b>Multi-Link</b> 5 lead, Multi-Link, Snap, 0.9m, AHA		10 pcs/box	
ELD-MQV5-90S-I	<b>Multi-Link</b> 5 lead, Multi-Link, Snap, 0.9m, IEC		10 pcs/box	
<b>3</b> ELD-MQV5-90G-A	<b>Multi-Link</b> 5 lead, Multi-Link, Grabber, 0.9m, AHA		10 pcs/box	
ELD-MQV5-90G-I	<b>Multi-Link</b> 5 lead, Multi-Link, Grabber, 0.9m, IEC		10 pcs/box	







## GE Healthcare

Part Number	Product Descriptions	OEM Cross Reference	Packaging	Price
1 ELD-TAF3-90G-A	<b>Apex Pro FH</b> 3 lead, Grabber, 0.9m, AHA		10 pcs/box	
ELD-TAF3-90G-I	<b>Apex Pro FH</b> 3 lead, Grabber, 0.9m, IEC		10 pcs/box	
2 ELD-TAF5-90G-A	<b>Apex Pro FH</b> 5 lead, Grabber, 0.9m, AHA		10 pcs/box	
ELD-TAF5-90G-I	<b>Apex Pro FH</b> 5 lead, Grabber, 0.9m, IEC		10 pcs/box	
ELD-TAF5-90S-A	<b>Apex Pro FH</b> 5 lead, Snap, 0.9m, AHA		10 pcs/box	
ELD-TAF5-90S-I	<b>Apex Pro FH</b> 5 lead, Snap, 0.9m, IEC		10 pcs/box	
ELD-TAF6-90S-A	<b>Apex Pro FH</b> 6 lead, Snap, 0.9m, AHA	33126	10 pcs/box	
ELD-TAF6-90S-I	<b>Apex Pro FH</b> 6 lead, Snap, 0.9m, IEC		10 pcs/box	
3 ELD-TAF6-90G-A	<b>Apex Pro FH</b> 6 lead, Grabber, 0.9m, AHA		10 pcs/box	
4 ELD-TAF6-90G-I	<b>Apex Pro FH</b> 6 lead, Grabber, 0.9m, IEC		10 pcs/box	



## Din Safety

Part Number	Product Descriptions	OEM Cross Reference	Packaging	Price
ELD-D3-60S-A	<b>Din Safety</b> 3 lead, Din Safety, Snap, 0.6m, AHA		10 pcs/box	
ELD-D3-60S-I	<b>Din Safety</b> 3 lead, Din Safety, Snap, 0.6m, IEC		10 pcs/box	
1 ELD-D3-60G-A	<b>Din Safety</b> 3 lead, Din Safety, Grabber, 0.6m, AHA		10 pcs/box	
ELD-D3-60G-I	<b>Din Safety</b> 3 lead, Din Safety, Grabber, 0.6m, IEC		10 pcs/box	
ELD-D5-60S-A	<b>Din Safety</b> 5 lead, Din Safety, Snap, 0.6m, AHA		10 pcs/box	
ELD-D5-60S-I	<b>Din Safety</b> 5 lead, Din Safety, Snap, 0.6m, IEC		10 pcs/box	
2 ELD-D5-60G-A	<b>Din Safety</b> 5 lead, Din Safety, Grabber, 0.6m, AHA		10 pcs/box	
ELD-D5-60G-I	<b>Din Safety</b> 5 lead, Din Safety, Grabber, 0.6m, IEC		10 pcs/box	
3 ELDE-D3-60-A	<b>Din Safety</b> 3 lead, Din Safety, with electrode, 0.6m, AHA		10 pcs/box	
4 ELDE-D3-60-I	<b>Din Safety</b> 3 lead, Din Safety, with electrode, 0.6m, IEC		10 pcs/box	
ELDE-D5-60-A	<b>Din Safety</b> 5 lead, Din Safety, with electrode, 0.6m, AHA		10 pcs/box	
ELDE-D5-60-I	<b>Din Safety</b> 5 lead, Din Safety, with electrode, 0.6m, IEC		10 pcs/box	
5 ELDER-D3-60-A	<b>Din Safety-Radiolucent</b> 3 lead, Din Safety, with electrode, 0.6m, AHA		10 pcs/box	
ELDER-D3-60-I	<b>Din Safety-Radiolucent</b> 3 lead, Din Safety, with electrode, 0.6m, IEC		10 pcs/box	
6 ELDER-D5-60-A	<b>Din Safety-Radiolucent</b> 5 lead, Din Safety, with electrode, 0.6m, AHA		10 pcs/box	
ELDER-D5-60-I	<b>Din Safety-Radiolucent</b> 5 lead, Din Safety, with electrode, 0.6m, IEC		10 pcs/box	



## Mindray-Datascope/Siemens

Part Number	Product Descriptions	OEM Cross Reference	Packaging	Price
ELD-DT3-90S-A	<b>Mindray-Datascope</b> 3 lead, Datascope, Snap, 0.9m, AHA	040-000748-01	10 pcs/box	
ELD-DT3-90S-I	<b>Mindray-Datascope</b> 3 lead, Datascope, Snap, 0.9m, IEC		10 pcs/box	
<b>1</b> ELD-DT3-90G-A	<b>Mindray-Datascope</b> 3 lead, Datascope, Grabber, 0.9m, AHA		10 pcs/box	
ELD-DT3-90G-I	<b>Mindray-Datascope</b> 3 lead, Datascope, Grabber, 0.9m, IEC		10 pcs/box	
ELD-DT5-90S-A	<b>Mindray-Datascope</b> 5 lead, Datascope, Snap, 0.9m, AHA	040-000746-01	10 pcs/box	
ELD-DT5-90S-I	<b>Mindray-Datascope</b> 5 lead, Datascope, Snap, 0.9m, IEC		10 pcs/box	
<b>2</b> ELD-DT5-90G-A	<b>Mindray-Datascope</b> 5 lead, Datascope, Grabber, 0.9m, AHA		10 pcs/box	
ELD-DT5-90G-I	<b>Mindray-Datascope</b> 5 lead, Datascope, Grabber, 0.9m, IEC		10 pcs/box	
ELD-SM5-90S-A	<b>Siemens</b> 5 lead, Dual Pin, Snap, 0.9m, AHA	33132	10 pcs/box	
ELD-SM5-90S-I	<b>Siemens</b> 5 lead, Dual Pin, Snap, 0.9m, IEC	33132E	10 pcs/box	
ELD-SM5-90G-A	<b>Siemens</b> 5 lead, Dual Pin, Grabber, 0.9m, AHA		10 pcs/box	
ELD-SM5-90G-I	<b>Siemens</b> 5 lead, Dual Pin, Grabber, 0.9m, IEC		10 pcs/box	





## Lead wires Adaptors

Part Number	Product Descriptions	OEM Cross Reference	Packaging	Price
1 ELA-TPH5-CV5	<b>Philips-Intellivue ECG Adaptor</b> 5 lead, used with Kendall/Covidien Leadwires	33526T,33526	1 pc/bag	
2 ELA-TMX5-PH5	<b>Philips-MX40 (ECG only) Adaptor</b> 5 lead, used with Philips-Intellivue Leadwires		1 pc/bag	
3 ELA-TMX5PH-PH5	<b>Philips-MX40 (ECG&amp; Philips Spo2) Adaptor</b> 5 lead, used with Philips-Intellivue Leadwires		1 pc/bag	
ELA-TMX5NE-PH5	<b>Philips-MX40 (ECG&amp; Nellcor Spo2) Adaptor</b> 5 lead, used with Philips-Intellivue Leadwires		1 pc/bag	
4 ELA-TMX5-MQ5	<b>Philips-MX40 (ECG only) Adaptor</b> 5 lead, used with GE Multi-LinkGE-Multi-Link Leadwires		1 pc/bag	
ELA-TMX5PH-MQ5	<b>Philips-MX40 (ECG&amp; Philips Spo2) Adaptor</b> 5 lead, used with GE Multi-LinkGE-Multi-Link Leadwires		1 pc/bag	
ELA-TMX5NE-MQ5	<b>Philips-MX40 (ECG&amp; Nellcor Spo2) Adaptor</b> 5 lead, used with GE Multi-LinkGE-Multi-Link Leadwires		1 pc/bag	
5 ELA-TMX5-CV5	<b>Philips-MX40 (ECG only) Adaptor</b> 5 lead, used with Kendall/Covidien Leadwires		1 pc/bag	
ELA-TMX5PH-CV5	<b>Philips-MX40 (ECG&amp; Philips Spo2) Adaptor</b> 5 lead, used with Kendall/Covidien Leadwires		1 pc/bag	
ELA-TMX5NE-CV5	<b>Philips-MX40 (ECG&amp; Nellcor Spo2) Adaptor</b> 5 lead, used with Kendall/Covidien Leadwires		1 pc/bag	



## Cardinal-Kendall


Part Number	Product Descriptions	OEM Cross Reference	Packaging	Price
ELD-CV3-90S-A	<b>Kendall</b> 3 lead, Snap, 0.9m, AHA	33103	10 pcs/box	
ELD-CV3-90S-I	<b>Kendall</b> 3 lead, Snap, 0.9m, IEC	33103E	10 pcs/box	
ELD-CV3-90G-A	<b>Kendall</b> 3 lead, Grabber, 0.9m, AHA		10 pcs/box	
ELD-CV3-90G-I	<b>Kendall</b> 3 lead, Grabber, 0.9m, IEC		10 pcs/box	
ELD-CV3-300S-A	<b>Kendall</b> 3 lead, Snap, 0.9m, AHA	33103	10 pcs/box	
ELD-CV3-300S-I	<b>Kendall</b> 3 lead, Snap, 0.9m, IEC	33103E	10 pcs/box	
ELD-CV3-300G-A	<b>Kendall</b> 3 lead, Grabber, 3.0m, AHA		10 pcs/box	
ELD-CV3-300G-I	<b>Kendall</b> 3 lead, Grabber, 3.0m, IEC		10 pcs/box	
<b>1</b> ELD-CV5-90S-A	<b>Kendall</b> 3 lead, Grabber, 0.9m, AHA	33105	10 pcs/box	
ELD-CV5-90S-I	<b>Kendall</b> 3 lead, Grabber, 0.9m, IEC	33105E	10 pcs/box	
<b>2</b> ELD-CV5-90G-A	<b>Kendall</b> 3 lead, Grabber, 0.9m, AHA		10 pcs/box	
ELD-CV5-90G-I	<b>Kendall</b> 3 lead, Grabber, 0.9m, IEC		10 pcs/box	
ELD-CV5-300S-A	<b>Kendall</b> 3 lead, Grabber, 3.0m, AHA	33105	10 pcs/box	
ELD-CV5-300S-I	<b>Kendall</b> 3 lead, Grabber, 3.0m, IEC	33105E	10 pcs/box	
ELD-CV5-300G-A	<b>Kendall</b> 3 lead, Grabber, 3.0m, AHA		10 pcs/box	
ELD-CV5-300G-I	<b>Kendall</b> 3 lead, Grabber, 3.0m, IEC		10 pcs/box	

**1**



## Cardinal-Kendall

Part Number	Product Descriptions	OEM Cross Reference	Packaging	Price
ELD-CVV5-90S-A	<b>Kendall</b> 3 lead, Grabber, 0.9m, AHA	33111	10 pcs/box	
ELD-CVV5-90S-I	<b>Kendall</b> 3 lead, Grabber, 0.9m, IEC	33111E	10 pcs/box	
① ETD-CV-CV5	<b>Kendall Trunk Cable</b> 5 lead, 2.2m, AHA/IEC	33145	10 pcs/box	
② ELDG-CV5-(22)90S-A	<b>Kendall Trunk Cable and Leadwire</b> 5 lead, Snap, 2.2m+0.9m, AHA	33135	10 pcs/box	
ELDG-CV5-(22)90S-I	<b>Kendall Trunk Cable and Leadwire</b> 5 lead, Snap, 2.2m+0.9m, IEC	33135E	10 pcs/box	
ELDG-CV5-(22)90G-A	<b>Kendall Trunk Cable and Leadwire</b> 5 lead, Grabber, 2.2m+0.9m, AHA		10 pcs/box	
ELDG-CV5-(22)90G-I	<b>Kendall Trunk Cable and Leadwire</b> 5 lead, Grabber, 2.2m+0.9m, IEC		10 pcs/box	



## Trunk Cables-Kendall/Covidien Leadwires

Part Number	Product Descriptions	OEM Cross Reference	Packaging	Price
1 ET-5005R(025)-CV-A	<b>Philips ECG Trunk Cable</b> 5 lead, 0.25m, used with Kendall/Covidien lead wires, AHA	33507	1 pc/bag	
ET-5005R(025)-CV-I	<b>Philips ECG Trunk Cable</b> 5 lead, 0.25m, used with Kendall/Covidien lead wires, IEC		1 pc/bag	
2 ET-5004(025)-CV-A	<b>GE Marquette ECG Trunk Cable</b> 5 lead, 0.25m, used with Kendall/Covidien lead wires, AHA	33506	1 pc/bag	
ET-5004(025)-CV-I	<b>GE Marquette ECG Trunk Cable</b> 5 lead, 0.25m, used with Kendall/Covidien lead wires, IEC		1 pc/bag	
ET-5016(025)-CV-A	<b>Mindray ECG Trunk Cable</b> 5 lead, 0.25m, used with Kendall/Covidien lead wires, AHA		1 pc/bag	
ET-5016(025)-CV-I	<b>Mindray ECG Trunk Cable</b> 5 lead, 0.25m, used with Kendall/Covidien lead wires, IEC		1 pc/bag	
ET-5016R(025)-CV-A	<b>Mindray ECG Trunk Cable(with Resister)</b> 5 lead, 0.25m, used with Kendall/Covidien lead wires, AHA		1 pc/bag	
ET-5016R(025)-CV-I	<b>Mindray ECG Trunk Cable(with Resister)</b> 5 lead, 0.25m, used with Kendall/Covidien lead wires, IEC		1 pc/bag	
ET-5038R(025)-CV-A	<b>Spacelabs ECG Trunk Cable</b> 5 lead, 0.25m, used with Kendall/Covidien lead wires, AHA		1 pc/bag	
ET-5038R(025)-CV-I	<b>Spacelabs ECG Trunk Cable</b> 5 lead, 0.25m, used with Kendall/Covidien lead wires, IEC		1 pc/bag	
3 ET-5007R(025)-CV-A	<b>Zoll ECG Trunk Cable</b> 5 lead, 0.25m, used with Kendall/Covidien lead wires, AHA	33509	1 pc/bag	
4 ET-5007R(025)-CV-I	<b>Zoll ECG Trunk Cable</b> 5 lead, 0.25m, used with Kendall/Covidien lead wires, IEC		1 pc/bag	

