

Single Patient Use NIBP Cuffs





Features




- Quality approved with FDA 510K clearance and CE mark
- Variety of connectors fitted with appropriate air hoses from a selection of brands
- Conformable Soft Fiber/Nonwoven BP Cuffs
- Wide range of patient sizes (from neonate to large adult)
- Single Patient Use eliminates cross-contamination issues




Every patient, every size,
every department
for all models of equipment

Neonate Single Patient Use Cuffs (Non-woven fabric)


Cuff Connector	Patient Size	Circumference	SKU	Packaging	Price
 <p>#1</p>	#1	3-6cm	BPD-11WS-W05P ¹	10 pcs/Box	
	#2	4-8cm	BPD-12WS-W05P	10 pcs/Box	
	#3	6-11cm	BPD-13WS-W05P	10 pcs/Box	
	#4	7-11cm	BPD-14WS-W05P	10 pcs/Box	
	#5	8-15cm	BPD-15WS-W05P	10 pcs/Box	
 <p>#2</p>	#1	3-6cm	BPD-11WS-W34	10 pcs/Box	
	#2	4-8cm	BPD-12WS-W34	10 pcs/Box	
	#3	6-11cm	BPD-13WS-W34	10 pcs/Box	
	#4	7-11cm	BPD-14WS-W34	10 pcs/Box	
	#5	8-15cm	BPD-15WS-W34	10 pcs/Box	
 <p>#3</p>	#1	3-6cm	BPD-11WS-W36 ²	10 pcs/Box	
	#2	4-8cm	BPD-12WS-W36	10 pcs/Box	
	#3	6-11cm	BPD-13WS-W36	10 pcs/Box	
	#4	7-11cm	BPD-14WS-W36	10 pcs/Box	
	#5	8-15cm	BPD-15WS-W36	10 pcs/Box	
 <p>#4</p>	#1	3-6cm	BPD-11WD-WD34	10 pcs/Box	
	#2	4-8cm	BPD-12WD-WD34	10 pcs/Box	
	#3	6-11cm	BPD-13WD-WD34	10 pcs/Box	
	#4	7-11cm	BPD-14WD-WD34	10 pcs/Box	
	#5	8-15cm	BPD-15WD-WD34	10 pcs/Box	
 <p>#5</p>	#1	3-6cm	BPD-11WD-WD05P ³	10 pcs/Box	
	#2	4-8cm	BPD-12WD-WD05P	10 pcs/Box	
	#3	6-11cm	BPD-13WD-WD05P	10 pcs/Box	
	#4	7-11cm	BPD-14WD-WD05P	10 pcs/Box	
	#5	8-15cm	BPD-15WD-WD05P	10 pcs/Box	




Standard Single Patient Use Cuffs (Non-woven fabric)

Cuff Connector	Patient Size	Circumference	SKU	Packaging	Price
 <p>Infant</p>	Infant	9-14.8cm	BPD-20WS-W01P	10 pcs/Box	
	Pediatric	13.8-21.5cm	BPD-30WS-W01P	10 pcs/Box	
	Small Adult	20.5-28.5cm	BPD-40WS-W01P	10 pcs/Box	
	Adult	27.5-36.5cm	BPD-50WS-W01P	10 pcs/Box	1
	Adult Long	27.5-49.5cm	BPD-60WS-W01P	10 pcs/Box	
	Large Adult	35.5-46cm	BPD-70WS-W01P	10 pcs/Box	
	Large Adult long	35.5-59cm	BPD-80WS-W01P	10 pcs/Box	
	Thigh	45-56.5cm	BPD-90WS-W01P	10 pcs/Box	
	 <p>Pediatric</p>	Infant	9-14.8cm	BPD-20WS-W08	10 pcs/Box
Pediatric		13.8-21.5cm	BPD-30WS-W08	10 pcs/Box	
Small Adult		20.5-28.5cm	BPD-40WS-W08	10 pcs/Box	
Adult		27.5-36.5cm	BPD-50WS-W08	10 pcs/Box	2
Adult Long		27.5-49.5cm	BPD-60WS-W08	10 pcs/Box	
Large Adult		35.5-46cm	BPD-70WS-W08	10 pcs/Box	
Large Adult long		35.5-59cm	BPD-80WS-W08	10 pcs/Box	
Thigh		45-56.5cm	BPD-90WS-W08	10 pcs/Box	
 <p>Small Adult</p>		Infant	9-14.8cm	BPD-20WS-W10	10 pcs/Box
	Pediatric	13.8-21.5cm	BPD-30WS-W10	10 pcs/Box	
	Small Adult	20.5-28.5cm	BPD-40WS-W10	10 pcs/Box	
	Adult	27.5-36.5cm	BPD-50WS-W10	10 pcs/Box	3
	Adult Long	27.5-49.5cm	BPD-60WS-W10	10 pcs/Box	
	Large Adult	35.5-46cm	BPD-70WS-W10	10 pcs/Box	
	Large Adult long	35.5-59cm	BPD-80WS-W10	10 pcs/Box	
	Thigh	45-56.5cm	BPD-90WS-W10	10 pcs/Box	


1





2




3




Standard Single Patient Use Cuffs (Non-woven fabric)

Cuff Connector	Patient Size	Circumference	SKU	Packaging	Price
 <p>Adult</p>	Infant	9-14.8cm	BPD-20WS-W11	10 pcs/Box	
	Pediatric	13.8-21.5cm	BPD-30WS-W11	10 pcs/Box	
	Small Adult	20.5-28.5cm	BPD-40WS-W11	10 pcs/Box	
	Adult	27.5-36.5cm	BPD-50WS-W11	10 pcs/Box	1
	Adult Long	27.5-49.5cm	BPD-60WS-W11	10 pcs/Box	
	Large Adult	35.5-46cm	BPD-70WS-W11	10 pcs/Box	
	Large Adult long	35.5-59cm	BPD-80WS-W11	10 pcs/Box	
	Thigh	45-56.5cm	BPD-90WS-W11	10 pcs/Box	
	 <p>Adult Long</p>	Infant	9-14.8cm	BPD-20WS-W37	10 pcs/Box
Pediatric		13.8-21.5cm	BPD-30WS-W37	10 pcs/Box	
Small Adult		20.5-28.5cm	BPD-40WS-W37	10 pcs/Box	
Adult		27.5-36.5cm	BPD-50WS-W37	10 pcs/Box	2
Adult Long		27.5-49.5cm	BPD-60WS-W37	10 pcs/Box	
Large Adult		35.5-46cm	BPD-70WS-W37	10 pcs/Box	
Large Adult long		35.5-59cm	BPD-80WS-W37	10 pcs/Box	
Thigh		45-56.5cm	BPD-90WS-W37	10 pcs/Box	




1




2




Standard Single Patient Use Cuffs (Non-woven fabric)

Cuff Connector	Patient Size	Circumference	SKU	Packaging	Price
 <p>Infant</p>	Infant	9-14.8cm	BPD-20WD-WD08	10 pcs/Box	
	Pediatric	13.8-21.5cm	BPD-30WD-WD08	10 pcs/Box	
	Small Adult	20.5-28.5cm	BPD-40WD-WD08	10 pcs/Box	
	Adult	27.5-36.5cm	BPD-50WD-WD08 ¹	10 pcs/Box	
	Adult Long	27.5-49.5cm	BPD-60WD-WD08	10 pcs/Box	
	Large Adult	35.5-46cm	BPD-70WD-WD08	10 pcs/Box	
	Large Adult long	35.5-59cm	BPD-80WD-WD08	10 pcs/Box	
	Thigh	45-56.5cm	BPD-90WD-WD08	10 pcs/Box	
	 <p>Pediatric</p>	Infant	9-14.8cm	BPD-20WD-W1011	10 pcs/Box
Pediatric		13.8-21.5cm	BPD-30WD-W1011	10 pcs/Box	
Small Adult		20.5-28.5cm	BPD-40WD-W1011	10 pcs/Box	
Adult		27.5-36.5cm	BPD-50WD-W1011 ²	10 pcs/Box	
Adult Long		27.5-49.5cm	BPD-60WD-W1011	10 pcs/Box	
Large Adult		35.5-46cm	BPD-70WD-W1011	10 pcs/Box	
Large Adult long		35.5-59cm	BPD-80WD-W1011	10 pcs/Box	
Thigh		45-56.5cm	BPD-90WD-W1011	10 pcs/Box	
 <p>Small Adult</p>		Infant	9-14.8cm	BPD-20WD-WD11	10 pcs/Box
	Pediatric	13.8-21.5cm	BPD-30WD-WD11	10 pcs/Box	
	Small Adult	20.5-28.5cm	BPD-40WD-WD11	10 pcs/Box	
	Adult	27.5-36.5cm	BPD-50WD-WD11 ³	10 pcs/Box	
	Adult Long	27.5-49.5cm	BPD-60WD-WD11	10 pcs/Box	
	Large Adult	35.5-46cm	BPD-70WD-WD11	10 pcs/Box	
	Large Adult long	35.5-59cm	BPD-80WD-WD11	10 pcs/Box	
	Thigh	45-56.5cm	BPD-90WD-WD11	10 pcs/Box	


1



2




3



4

Standard Single Patient Use Cuffs (Non-woven fabric)

Cuff Connector	Patient Size	Circumference	SKU	Packaging	Price
 <p>Adult</p>	Infant	9-14.8cm	BPD-20WD-W38	10 pcs/Box	
	Pediatric	13.8-21.5cm	BPD-30WD-W38	10 pcs/Box	
	Small Adult	20.5-28.5cm	BPD-40WD-W38	10 pcs/Box	
	Adult	27.5-36.5cm	BPD-50WD-W38	10 pcs/Box	¹
	Adult Long	27.5-49.5cm	BPD-60WD-W38	10 pcs/Box	
	Large Adult	35.5-46cm	BPD-70WD-W38	10 pcs/Box	
	Large Adult long	35.5-59cm	BPD-80WD-W38	10 pcs/Box	
	Thigh	45-56.5cm	BPD-90WD-W38	10 pcs/Box	

¹

