

# Incubator/Warmer Temperature Probes



## Features

- *Specially design for neonate skin temperature measurement*
- *Engineered to function equivalent to OEMs*
- *Quality approved with CE Certificate and FDA 510K*
- *Accuracy:  $\pm 0.1$  C*
- *Affordable prices*
- *12 months warranty*




**New replacements for specific models by these major brands**

- **Air Shields/Hill Rom/Drager**
- **Atom**
- **David**
- **Dison**
- **Ohmeda**


## Incubator/Warmer Temperature Probes

Part Number	Product Descriptions	OEM Cross Reference	Packaging	Price
1 TS-DGA-PS15	<b>Air Shields/Hill-Rom/Draeger</b> Reusable Skin Temperature Probe, For Isolette C2000 & 8000, Versalet 7700	MU12533, 83-102-41	1 pc/bag	
TSD-DGA-GS15	<b>Air Shields/Hill-Rom/Draeger</b> Single-use Skin Temperature Probe with reflective cover, For Isolette C2000 & 8000, Versalet 7700	MU12551, 83-104-15	25 pcs/box	
2 TS-DGC-PS15	<b>Air Shields/Hill-Rom/Draeger</b> Reusable Skin Temperature Probe, For Resuscitaire (RW82-1)	MU11225, 81-300-05	1 pc/bag	
TSD-DGC-GS15	<b>Air Shields/Hill-Rom/Draeger</b> Single-use Skin Temperature Probe with reflective cover, For Resuscitaire(RW82-1)	MU11228, 81-300-08	25 pcs/box	
3 TS-AS-PS15	<b>Air Shields/Hill-Rom/Draeger</b> Reusable Skin Temperature Probe, For IICS(PM78-1), Isolette C-100, C-100QT, C-450QT & 550QT	MU06951, 68-209-70	1 pc/bag	
TSD-AS-GS15	<b>Air Shields/Hill-Rom/Draeger</b> Single-use Skin Temperature Probe with reflective cover, For IICS(PM78-1), Isolette C-100, C-100QT, C-450QT & 550QT	MU06933, 68-209-20	25 pcs/box	
4 TS-DGB-PS15	<b>Air Shields/Hill-Rom/Draeger</b> Reusable Skin Temperature Probe, For Babytherm8000 WB/8004/8010/4200	MX11000	1 pc/bag	
TSD-DGB-GS15	<b>Air Shields/Hill-Rom/Draeger</b> Single-use Skin Temperature Probe with reflective cover, For Babytherm8000 WB/8004/8010/4200	MX11000	25 pcs/box	



## Incubator/Warmer Temperature Probes

Part Number	Product Descriptions	OEM Cross Reference	Packaging	Price
1 TS-OHA-PS15	<b>Ohmeda</b> Reusable Skin Temperature Probe, For Care Plus, GC, IC & IWS	6600-0628-700, 0208-0697-700	1 pc/bag	
TSD-OHA-GS15	<b>Ohmeda</b> Single-use Skin Temperature Probe with reflective cover, For Care Plus, GC, IC & IWS	6600-0208-700	25 pcs/box	
2 TS-OH-PS15	<b>Ohmeda</b> Reusable Skin Temperature Probe, For Giraffe Omnibed	6600-0875-700	1 pc/bag	
TSD-OH-GS15	<b>Ohmeda</b> Single-use Skin Temperature Probe with reflective cover, For Giraffe Omnibed	6600-0873-700, 6600-0874-700	25 pcs/box	
TS-AT-PS15	<b>Atom</b> Reusable Skin Temperature Probe, For Model 100/101/103, V-2100/2100G/2200/550/505ST/808	NF-024(yellow)	1 pc/bag	
3 TSD-AT-GS15	<b>Atom</b> Single-use Skin Temperature Probe with reflective cover, For Model 100/101/103, V-2100/2100G/2200/550/505ST/808	NF-023(yellow), 60881(white)	25 pcs/box	
4 TS-FA-PS15	<b>FANEM</b> Reusable skin temp. probe		1 pc/bag	



## Incubator/Warmer Temperature Probes

Part Number	Product Descriptions	OEM Cross Reference	Packaging	Price
1 TS-DVA-PS15	<b>David</b> Reusable Skin Temperature Probe, For HKN-93A, Resuscitairre Infant, Radiant warmer		1 pc/bag	
2 TS-DVB-PS15	<b>David</b> Reusable Skin Temperature Probe, for TI-2000 transport incubator		1 pc/bag	
3 TS-DVC-PS15	<b>David</b> Reusable Skin Temperature Probe, For Incubators and Infant warmers		1 pc/bag	
4 TS-DSA-PS15	<b>Dison</b> Reusable Skin Temperature Probe	810040007	1 pc/bag	
5 TS-DSB-PS15	<b>Dison</b> Reusable Skin Temperature Probe, For all requied as second temperature probe	810040002	1 pc/bag	
6 TS-BE01-PS15	<b>Being</b> Reusable Skin Temperature Probe, For Redel 2 pins connector		1 pc/bag	
7 TS-BE02-PS15	<b>Being</b> Reusable Skin Temperature Probe, Fo Redel 4 pins connector		1 pc/bag	

