

## Humidification Chamber

The reusable brand-new replacement parts for Ventilators

### Order Information



Adult Type



Pediatric/Infant Type

Part Number	Patient Sizes	Gas in Connector	Gas out Connector	Price (USD)
EHC-A01	Adult	ID 22mm OD 27mm	ID 15mm OD 22mm	\$50.00
EHC-P01	Pediatric/Infant	ID 15mm OD 22mm	ID 15mm OD 22mm	\$55.00

### Specifications

Cap Material	Silicone	Sterilize	Autoclavable 134 °C
Bottle Material	PC medical grade	Bottom Diameter	124.5 mm
Height(Adult)	115 mm	Height ( Pediatric/Infant )	70 mm
Bottle Color	Clear	Packaging Type	Carton box
Packaging Unit	1pcs/box	Base Material	Aluminum
O Ring Material	Silicone	Patient Size	4kg above
Warranty	12 months		

